

ASEPTI ACTIVE WIPES

Hospital grade disinfectant

ARTG 182376
ARTG 182377

Asepti Active Wipes offer hospital grade disinfectant that meets the TGA requirement for a single step clean and disinfect suitable to be used in “dirty conditions”.



DISINFECTANT CLEANER

Asepti Active Wipes contain three actives including two latest technology QAC's effective against traditionally QAC resistant organisms.

Asepti Active Wipes provide broad spectrum bactericidal efficacy against both gram positive and gram negative organisms including Methicillin-Resistant Staphylococcus Aureus (MRSA) and Vancomycin-Resistant Enterococcus (VRE)

Asepti Active Wipes were evaluated in accordance with the TGA test for Hospital Grade Disinfectants under "dirty conditions", thereby approved for use as a "single step" clean and disinfect.

Asepti Active Wipes' ready to use formats are simple and easy to use, they ensure the correct concentration of chemistry is available in each application and reduce the risk of cross contamination.

Asepti Active Wipes are conveniently packed and easy to use, no lint, no foam.

Asepti Active Wipes are compatible with **Asepti Active Liquid**.

A range of dispensers to suit the flat pack and pail format are also available encouraging accessibility and promoting usage.

DIRECTIONS FOR USE

Asepti Active Wipes are safe for use on most surfaces. They are not suitable for therapeutic devices or as a baby wipe.

Asepti Active Wipes are single use only, supporting infection control protocols.

SAFETY INFORMATION

Do not mix with any other products. For further information and safety instructions see product label and Safety Data Sheet.

Product Code	Product Description	Pack Size
7101608	Asepti Active Wipes	12 x 100 wipes
99960090	Flat Pack Holder	Each
7100096	Asepti Active Wipes	3 x 325 wipes
9996-1208	325 Wipe Wall Bracket	Each
7100752	Asepti Active Liquid	12 x 500mL
401512	Wall and bed bracket	Each

PLEASE SPEAK TO YOUR HEALTHCARE MANAGER OR CALL **1800 022 002** FOR FURTHER INFORMATION